


OPEN SURGERY



Wundspreizer

Self-Retaining Retractors

Separadores autoestáticos









Ecarteurs autostatiques

Divaricatori autostatici

070

Wundspreizer
Self-Retaining Retractors

Separadores autoestáticos
Ecarteurs autostatiques
Divaricatori autostatici

Art.Nr.	Beschreibung Description	Size	Spez.
GELPI			
070-202-09	9,0 cm; 3 1/2"		
070-202-11	11,0 cm; 4 3/8"		
070-202-14	13,5 cm; 5 3/8"		
GELPI			
070-202-19	19,0 cm; 7 4/8"		
DE BAKEY			
070-240-10	10,0 cm; 3 7/8"		
070-240-14	14,0 cm; 5 4/8"		
GELPI-LOKTITE			
070-250-19	19,0 cm; 7 4/8"		


Wundspreizer
Self-Retaining Retractors

Separadores autoestáticos


Ecarteurs autostatiques

Divaricatori autostatici

Art.Nr.	Beschreibung Description	Size	Spez.
070-300-00	HENLEY	16,5 cm; 6 4/8"	SET



scharf | sharp | agudo | pointu | acuto



1/1




- bestehend aus:**
- 070-300-16 1 Sperrer, scharf
 - 070-300-19 1 gezahnte Mittelvalve 19 x 16 mm
 - 070-300-25 1 gezahnte Mittelvalve 25 x 16 mm
 - 070-300-32 1 gezahnte Mittelvalve 30 x 16 mm

- consisting of:**
- 070-300-16 1 Spreader, sharp
 - 070-300-19 1 Central valve, serrated, 19 x 16 mm
 - 070-300-25 1 Central valve, serrated, 25 x 16 mm
 - 070-300-32 1 Central valve, serrated, 30 x 16 mm

- compuesto de:**
- 070-300-16 1 Separador, agudo
 - 070-300-19 1 Valva central, dentada, 19 x 16 mm
 - 070-300-25 1 Valva central, dentada, 25 x 16 mm
 - 070-300-32 1 Valva central, dentada, 30 x 16 mm

- composé de:**
- 070-300-16 1 Ecarteur, griffes pointues
 - 070-300-19 1 Valve centrale, dentelée, 19 x 16 mm
 - 070-300-25 1 Valve centrale, dentelée, 25 x 16 mm
 - 070-300-32 1 Valve centrale, dentelée, 30 x 16 mm

- composto da:**
- 070-300-16 1 Divaricatore, acuto
 - 070-300-19 1 Valva centrale, dentelata, 19 x 16 mm
 - 070-300-25 1 Valva centrale, dentelata, 25 x 16 mm
 - 070-300-32 1 Valva centrale, dentelata, 30 x 16 mm



Wundspreizer
Self-Retaining Retractors

Separadores autoestáticos

Ecarteurs autostatiques


Divaricatori autostatici

TONTARRA

Art.Nr.	Beschreibung Description	Size	Spez.
070-301-00	HENLEY	16,5 cm; 6 4/8"	SET



stumpf | blunt | romo | mousse | smusso



1/1




- bestehend aus:**
- 070-301-16 1 Sperrer, stumpf
 - 070-300-19 1 gezahnte Mittelvalve 19 x 16 mm
 - 070-300-25 1 gezahnte Mittelvalve 25 x 16 mm
 - 070-300-32 1 gezahnte Mittelvalve 30 x 16 mm

- consisting of:**
- 070-301-16 1 Spreader, blunt
 - 070-300-19 1 Central valve, serrated, 19 x 16 mm
 - 070-300-25 1 Central valve, serrated, 25 x 16 mm
 - 070-300-32 1 Central valve, serrated, 30 x 16 mm

- compuesto de:**
- 070-301-16 1 Separador, romo
 - 070-300-19 1 Valva central, dentada, 19 x 16 mm
 - 070-300-25 1 Valva central, dentada, 25 x 16 mm
 - 070-300-32 1 Valva central, dentada, 30 x 16 mm

- composé de:**
- 070-301-16 1 Ecarteur, mousse
 - 070-300-19 1 Valve centrale, dentelée, 19 x 16 mm
 - 070-300-25 1 Valve centrale, dentelée, 25 x 16 mm
 - 070-300-32 1 Valve centrale, dentelée, 30 x 16 mm



- composto da:**
- 070-301-16 1 Divaricatore, smusso
 - 070-300-19 1 Valva centrale, dentelata, 19 x 16 mm
 - 070-300-25 1 Valva centrale, dentelata, 25 x 16 mm
 - 070-300-32 1 Valva centrale, dentelata, 30 x 16 mm



Wundspreizer
Self-Retaining Retractors

Separadores autoestáticos
Ecarteurs autostatiques
Divaricatori autostatici

Art.Nr.	Beschreibung Description	Size	Spez.
070-320-17	MAYO-ADAMS mit Mittelvalve with central blade con valva central avec valve centrale con valva centrale	17,0 cm; 6 6/8"	




Wundspreizer
Self-Retaining Retractors

Separadores autoestáticos
Ecarteurs autostatiques
Divaricatori autostatici


Art.Nr. Beschreibung | Description Size Spez.

MEYERDING

20,0 cm; 7 7/8"



070-482-04




65 x 25 mm




65 x 25 mm

070-482-02




55 x 25 mm




55 x 25 mm

070-482-01



32 x 25 mm



32 x 25 mm

Art.Nr. Beschreibung | Description Size Spez.


FINSEN



070-610-05	scharf sharp agudo pointu acuto	2 x 3		5,0 cm; 2"
070-611-05	stumpf blunt romo mousse smusso	2 x 3		5,0 cm; 2"
070-610-07	scharf sharp agudo pointu acuto	3 x 4		7,0 cm; 2 6/8"
070-611-07	stumpf blunt romo mousse smusso	3 x 4		7,0 cm; 2 6/8"

1/1

ALM



070-620-07	scharf sharp agudo pointu acuto	4 x 4		7,0 cm; 2 6/8"
070-622-07	halbscharf semi-sharp semi-agudo semi-pointu semi-acuto	4 x 4		7,0 cm; 2 6/8"
070-621-07	stumpf blunt romo mousse smusso	4 x 4		7,0 cm; 2 6/8"
070-620-10	scharf sharp agudo pointu acuto	4 x 4		10,0 cm; 3 7/8"
070-622-10	halbscharf semi-sharp semi-agudo semi-pointu semi-acuto	4 x 4		10,0 cm; 3 7/8"
070-621-10	stumpf blunt romo mousse smusso	4 x 4		10,0 cm; 3 7/8"

1/1

Wundspreizer
Self-Retaining Retractors

Separadores autoestáticos
Ecarteurs autostatiques
Divaricatori autostatici

Art.Nr.	Beschreibung Description	Size	Spez.
---------	----------------------------	------	-------


JANSEN

10,0 cm; 3 7/8"




070-630-10 scharf | sharp
agudo | pointu | acuto

3 x 3



070-631-10 stumpf | blunt
romo | mousse | smusso

3 x 3



1/1


SCALP-CONTOUR

14,0 cm; 5 4/8"




070-660-14 scharf | sharp
agudo | pointu | acuto

3 x 4



070-661-14 stumpf | blunt
romo | mousse | smusso

3 x 4





1/1

070

Wundspreizer
Self-Retaining Retractors

Separadores autoestáticos
 Ecarteurs autostatiques
 Divaricatori autostatici

Art.Nr.	Beschreibung Description	Size	Spez.
	WEITLANER-LOKTITE	10,5 cm; 4 1/8"	



070-670-10

scharf | sharp
 agudo | pointu | acuto


2 x 3



070-671-10

stumpf | blunt
 romo | mousse | smusso

2 x 3



1/1

070

Wundspreizer
Self-Retaining Retractors

Separadores autoestáticos
Ecarteurs autostatiques
Divaricatori autostatici


Art.Nr. Beschreibung | Description Size Spez.

WEITLANER




070-670-13 scharf | sharp 3 x 4 13,0 cm; 5 1/8"

agudo | pointu | acuto




070-671-13 stumpf | blunt 3 x 4 13,0 cm; 5 1/8"

romo | mousse | smusso




070-670-16 scharf | sharp 3 x 4 16,5 cm; 6 4/8"

agudo | pointu | acuto



070-671-16 stumpf | blunt 3 x 4 16,5 cm; 6 4/8"

romo | mousse | smusso




1/1

070-674-15 **MOLLISON** 15,5 cm; 6 1/8"

scharf | sharp | agudo | pointu | acuto



4 x 4




1/1

Wundspreizer
Self-Retaining Retractors

Separadores autoestáticos
Ecarteurs autostatiques
Divaricatori autostatici

Art.Nr. Beschreibung | Description Size Spez.


WEITLANER



070-670-20

scharf | sharp
agudo | pointu | acuto

3 x 4




20,0 cm; 7 7/8"

070-671-20

stumpf | blunt
romo | mousse | smusso

3 x 4




20,0 cm; 7 7/8"

070-670-24

scharf | sharp
agudo | pointu | acuto

3 x 4




24,0 cm; 9 4/8"

070-671-24

stumpf | blunt
romo | mousse | smusso

3 x 4



24,0 cm; 9 4/8"

Wundspreizer
Self-Retaining Retractors

Separadores autoestáticos
Ecarteurs autostatiques
Divaricatori autostatici

Art.Nr. Beschreibung | Description Size Spez.



PLESTER

11,0 cm; 4 3/8"



scharf | sharp | agudo | pointu | acuto



070-672-11





070-673-11



1/1

PLESTER

13,0 cm; 5 1/8"



070-830-13

scharf | sharp
agudo | pointu | acuto

2 x 2



070-831-13

stumpf | blunt
romo | mousse | smusso

2 x 2




1/1

Wundspreizer
Self-Retaining Retractors


Separadores autoestáticos
 Ecarteurs autostatiques
 Divaricatori autostatici

Art.Nr.	Beschreibung Description	Size	Spez.
	ANDERSON-ADSON	20,0 cm; 7 7/8"	




070-678-20 **scharf | sharp**
 agudo | pointu | acuto

4 x 4



070-679-20 **stumpf | blunt**
 romo | mousse | smusso

4 x 4



1/1


070

Wundspreizer
Self-Retaining Retractors


Separadores autoestáticos
Ecarteurs autostatiques
Divaricatori autostatici

Art.Nr. Beschreibung | Description Size Spez.


ADSON




070-682-26 scharf sharp agudo pointu acuto 3 x 4 26,0 cm; 10 2/8"



070-683-26 stumpf blunt romo mousse smusso 3 x 4 26,0 cm; 10 2/8"




070-682-29 scharf sharp agudo pointu acuto 6 x 6 29,0 cm; 11 3/8"



à

070-683-29 stumpf blunt romo mousse smusso 6 x 6 29,0 cm; 11 3/8"



Wundspreizer
Self-Retaining Retractors


Separadores autoestáticos
Ecarteurs autostatiques
Divaricatori autostatici

Art.Nr.	Beschreibung Description	Size	Spez.
	ADSON	28,5 cm; 11 2/8"	
	scharf sharp agudo pointu acuto		





070

070-686-01



6 x 6

lange Zähne rechts
long teeth on right side
dientes largos a la derecha
griffes longues à droite
denti lunghe a destra





070-686-02



6 x 6

lange Zähne links
long teeth on left side
dientes largos a la izquierda
griffes longues à gauche
denti lunghe a sinistra




Wundspreizer
Self-Retaining Retractors

Separadores autoestáticos
Ecarteurs autostatiques
Divaricatori autostatici

Art.Nr.	Beschreibung Description	Size	Spez.
070-690-31	BECKMANN	31,0 cm; 12 2/8"	


scharf | sharp | agudo | pointu | acuto



1/1

070-691-32	ADSON	32,0 cm; 12 5/8"	
------------	--------------	------------------	--

stumpf | blunt | romo | mousse | smusso



1/1

Wundspreizer
Self-Retaining Retractors

Separadores autoestáticos

Ecarteurs autostatiques

Divaricatori autostatici

Art.Nr.	Beschreibung Description	Size	Spez.
	WILLIAMS	17,0 cm; 6 6/8"	




Blattgröße Blade size Tamaño hoja Taille de valve Misura lame	links left izquierdo gauche sinistro	rechts right derecho droite destro
40 x 10 mm	070-700-17	070-710-17
50 x 10 mm	070-701-17	070-711-17
60 x 10 mm	070-702-17	070-712-17
70 x 10 mm	070-703-17	070-713-17
50 x 20 mm	070-704-17	070-714-17
60 x 20 mm	070-705-17	070-715-17
70 x 20 mm	070-706-17	070-716-17


Wundspreizer
Self-Retaining Retractors

Separadores autoestáticos
Ecarteurs autostatiques
Divaricatori autostatici


Art.Nr.	Beschreibung Description	Size	Spez.
	WULLSTEIN		13,0 cm; 5 1/8"



070-810-13 **scharf | sharp** **3 x 3**
agudo | pointu | acuto




070-811-13 **stumpf | blunt** **3 x 3**
romo | mousse | smusso




1/1

070-840-13	MILLIGAN	13,0 cm; 5 1/8"
	stumpf blunt romo mousse smusso	









3 x 3



1/1

Wundspreizer
Self-Retaining Retractors

Separadores autoestáticos
Ecarteurs autostatiques
Divaricatori autostatici


Art.Nr.	Beschreibung Description		Size	Spez.
ADSON				
				
070-860-13	scharf sharp agudo pointu acuto	3 x 4		13,5 cm; 5 3/8"
070-860-16	scharf sharp agudo pointu acuto	3 x 4		16,0 cm; 6 2/8"
070-861-13	stumpf blunt romo mousse smusso	3 x 4		13,5 cm; 5 3/8"
070-861-16	stumpf blunt romo mousse smusso	3 x 4		16,0 cm; 6 2/8"

Wundspreizer
Self-Retaining Retractors

Separadores autoestáticos
Ecarteurs autostatiques
Divaricatori autostatici

Art.Nr. Beschreibung | Description Size Spez.


ADSON



070-860-20

scharf | sharp
agudo | pointu | acuto

3 x 4




20,0 cm, 7 7/8"

070-860-26

scharf | sharp
agudo | pointu | acuto

3 x 4




26,0 cm; 10 2/8"

070-861-20

stumpf | blunt
romo | mousse | smusso

3 x 4




20,0 cm; 7 7/8"

070-861-26

stumpf | blunt
romo | mousse | smusso

3 x 4



26,0 cm; 10 2/8"


070

Wundspreizer
Self-Retaining Retractors


Separadores autoestáticos
Ecarteurs autostatiques
Divaricatori autostatici

Art.Nr.	Beschreibung Description	Size	Spez.
070-862-32	BECKMANN-EATON	32,0 cm; 12 5/8"	

scharf | sharp | agudo | pointu | acuto




7 x 7



1/1

070-850-13	WEITLANER-BABY	13,5 cm; 5 3/8"	
------------	-----------------------	-----------------	--



1/1

Wundspreizer
Self-Retaining Retractors

Separadores autoestáticos
Ecarteurs autostatiques
Divaricatori autostatici

Art.Nr. Beschreibung | Description Size Spez.


BECKMANN-ADSON

31,0 cm; 12 2/8"


scharf | sharp | agudo | pointu | acuto




070-870-01 25 x 25 mm 4 x 4




070-870-02 25 x 31 mm 4 x 4



070-870-03 25 x 40 mm 4 x 4



070-870-04 25 x 45 mm 4 x 4



Wundspreizer
Self-Retaining Retractors


Separadores autoestáticos
Ecarteurs autostatiques
Divaricatori autostatici

Art.Nr.	Beschreibung Description	Size	Spez.
	BECKMANN-ADSON		31,0 cm; 12 2/8"


stumpf | blunt | romo | mousse | smusso




070-871-01 25 x 25 mm 4 x 4



070-871-02 25 x 31 mm 4 x 4



070-871-03 25 x 40 mm 4 x 4



070-871-04 25 x 45 mm 4 x 4

